

FRAXPEN

BEAUTY & FACE

Technical Specifications:

Dimension: 177.5mm x 34.88mm x 34.88mm
Cradle: 80.65mm x 80.65mm x 68.21mm
Power supply: Battery Lithium polymer, 500mAh
Weight: 530g
Duration: Recharging time: 3 hrs Continuous using time: 3 hrs
High speed mode: 700 cycles/min Low speed mode: 412 cycles/min



ABALASE
ENERGY INNOVATIONS

22 Itis str., Athens
Greece 16344

t +30 2109765002-6
f +30 2109765001

abalase@otenet.gr
www.abalase.gr

ABALASE
ENERGY INNOVATIONS



FRAXPEN
NANO ENERGY

CE



Treatment Principle

Fraxpen increases e-conductivity fibroblast growth factors, transforming growth factors, platelet derived growth factors etc, obtaining new collagen and elastin production (collagen III and collagen I). Vibration and microstimulation helps proliferation of fibroblasts. Fast healing in less than an hour with minimal injury.

Technology

Fraxpen uses electronic automated vibro heads tips that move perpendicularly from the surface to the skin. Vibration wave on the tip is 700 times per minute, maximum creating 6.300 micro channels through the skin.

Easy to use, more time-efficient and less painful than the micro roller heads. The tip has nano vibro needles which are efficient to use in sensitive areas such as around nose or eyes.

Indications

Face

Skin Rejuvenation
Enhancement of Skin Elasticity
Acne Scar Reduction
Whitening/ Brightening

Eyes

Wrinkle Reduction
Skin toning
Tightening/ Lifting
Eyelid area treatment

Special Features

- Automatic and rechargeable battery. Easy and safe to use because of no needles.
- The depth of nano needle trauma is adjustable. You only have to change the guide if you need shorter or longer vibration for different areas or purpose of the treatment.
- 5 different speed modes provide both safety and convenience. The Low Speed Mode allows slow and careful treatment without the risk of excessive irritation to the skin. If the procedure needs to be more time-efficient, high Speed Modes create adequate number of micro channels even with the faster move.
- All Modes create channels with proper spaces between them, without micro needling. This is the reason that Fraxpen is able to cover the same area in shorter time.
- It is less painful. The vibro tip is touching your skin at a perpendicularly, way. Pain is minimized on no pain at all.
- It is very easy to treat small areas such as around nose, eyes, or mouth. The hand-piece allows to treat areas in whichever angle you wish with simultaneous emission of mesotherapy solutions.