

SYNERGAN

When Skin Health meets Technology


Technical Specs

Work Mode:	Hydporation, CO ₂ Therapy, MultiPolar RF, Ultrasound
Hydporation Parameters:	Vacuum Scope 25 kPa-90kPa Bubble Size: 10 μm
CO ₂ Therapy Parameters:	Energy Level: 10-100; Increment: 10
MultiPolar RF Parameters:	Electrode Quantity: 4 Energy Level: 10-100, Increment: 10 Frequency: MHz
Ultrasound Parameters:	Ultrasound Energy Level: 10-100, Increment: 10, Ultrasound Frequency: 1MHz Vibration Energy Level: 10-100, Increment: 10, Increment: 10
Rated Input Power:	50VA
Net Weight:	20Kg
Machine Size:	L 500mm x W 460mm x H 1050mm



Professional Facial and Body Therapies.
Anti-Aging Protocols.
4 Technologies in one device.

- Hydor • CO₂
- RF4 • US4



ABALASE

22, Itis Str., 16344, Athens, Greece
T 210 9765002
F 210 9765001
e abalase@otenet.gr



SYNERGAN

When Skin Health meets Technology

- Deep hydrate and retains the Skin's Moisture
- Increases the thickness of the skin
- Restores the volume of the skin
- Renews and revives a youthful glow
- Reduces pigmentation and discoloration of the skin
- Improves cellular turnover
- Eliminates the appearance of fine lines and wrinkles
- Reduces the static rhytides
- Increase skin's firmness and elasticity
- Refines skin's texture
- Maintain an even skin tone
- Soothes, calms and provides anti-oxidant protection during summer time



CO₂

It is based on the Bohr Effect, a physiological phenomenon. High Concentrations of CO₂ trigger cellular response for oxygenation. The skin is naturally regenerated. It creates an optimal environment for the infusion of essential nutrients, both during and after treatments.



RF4

The RF4 technology is a new, revolutionary skin rejuvenating system, giving instant results. This controlled RF energy delivers instant visible results because of the flow of electrons. The treatment targets on regeneration of collagen fibers in the dermis, resulting in firming and contouring the oval structure of the face. Ideal for the elimination of fine lines and rhytides.



HYDOR

Hydor, a more advanced technique generates a power to form a continuous liquid flow on the epidermis. It cleanses impurities and dead skin cells, thus it is clarifying the complexion by promoting cellular turnover in the most safe and painless way. Hydor, gives the skin healthy glow, and prepares the skin for the next stages of innovative therapies.



US4

The stimulating vibration increases the blood circulation in the capillaries and activates the metabolic functions. The US4 technology provides deeper penetration and absorption of essential minerals and vitamins. The water molecules will be easily broken into smaller ones carrying all nutrients deep into skin layers.