

AKOUΩ

ACOUSTIC WAVE THERAPY

Indications:

- Cellulite
- Skin tightening
- Skin firmness
- Stretch marks
- Lymphatic drainage
- Body shaping

Technical Specifications:

- Energy/ Intensity: 60-185mJ (1-5bars)
- Working frequency: 1-22Hz
- Work pattern: 4 continuous capture modes (shots 4/8/12 pulses)
- Handle description: anodized aluminum shell with fan cooling
- Power consumption: 100/240VAC 50/60 Hz 300VA
- Wire fuse specification: 250V/10A
- Screen: Color touch screen
- Dimensions: 290x240x130mm (length/width/height)
- Product weight: 2.07kg (do not contain therapeutic handle & foot switch)



AKOUΩ

ACOUSTIC WAVE THERAPY

Quick & Easy

AKOUΩ


ACOUSTIC WAVE THERAPY

The best and fastest for:

Cellulite
Skin tightening
Stretch marks
Lymphatic drainage

Combinations

- Cryolipolysis
- RF
- Ultrasound



Technique principle:

Acoustic Waves penetrates into different layers of the skin repairing the cellulite-causing connective tissue. Break down the fibrous walls releasing the clogged fat cells. Increase blood flow of treated areas and helps the lymphatic system for better detoxification. Once the circulation is stimulated, the broken down fat cells and toxins are directed to the liver naturally and metabolised. As a result, you will get smoother, firmer and more youthful skin after only a few sessions.

Mechanism of Action

- Improves blood circulation
- Improves lymphatic circulation
- Improves intracellular liquid
- Minimizes swelling and fluid retention
- Activates cellular regeneration
- Improves skin's complexion
- Strengthens connective tissue
- Reduces the adipocyte volume (lipolysis)
- Liposculpture

Advantages

- Safe
- Transportable
- Non invasive
- Pain free
- Quick & easy
- No down time
- Immediate results

